# BD Oncomark CD14 FITC/CD64 PE

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7/2010

23-13735-01

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1. Olweus J, Terstappen LWMM, Thompson PA, Lund-Johansen F. Expression and function of receptors for stem cell factor and erythropoietin during lineage commitment of human hematopoietic progenitor cells. *Blood*. 1996;88:1594-1607.
2. Krasinskas AM, Wasik MA, Kamoun M, Schretzenmair R, Moore J, Salhany K. The usefulness of CD64, other monocyte-associated antigens, and CD45 gating in the subclassification of acute myeloid leukemias with monocytic differentiation. *Am J Clin Pathol*. 1998;110:797-805.
3. Dimitriu-Bona A, Burmester GR, Waters SJ, Winchester RJ. Human mononuclear phagocyte differentiation antigens. I. Patterns of antigenic expression on the surface of human monocytes and macrophages defined by monoclonal antibodies. *J Immunol*. 1983;130:145-152.
4. Herrmann F, Komischke B, Odenwald E, Ludwig WD. Use of monoclonal antibodies as a diagnostic tool in human leukemia. I. Acute myeloid leukemia and acute phase of chronic myeloid leukemia. *Blut*. 1983;47:157-163.
5. Rothe G, Schmitz G. Consensus protocol for the flow cytometric immunophenotyping of hematopoietic malignancies. *Leukemia*. 1996;10:877-895.
6. Stelzer GT, Marti G, Hurley A, McCoy P, Jr., Lovett EJ, Schwartz A. US-Canadian consensus recommendations on the immunophenotypic analysis of hematologic neoplasia by flow cytometry: standardization and validation of laboratory procedures. *Cytometry*. 1997;30:214-230.
7. *Protection of Laboratory Workers from Occupationally Acquired Infections—Second Edition; Approved Guideline*. Villanova, PA: National Committee for Clinical Laboratory Standards; 2001. NCCLS document M29-A2.
8. Centers for Disease Control. Update: universal precautions for prevention of transmission of human immunodeficiency virus, hepatitis B virus, and other bloodborne pathogens in health-care settings. *MMWR*. 1988;37:377-388.
9. Jackson AL, Warner NL. Preparation, staining, and analysis by flow cytometry of peripheral blood leukocytes. In: Rose NR, Friedman H, Fahey JL, eds. *Manual of Clinical Laboratory Immunology*. 3rd ed. Washington, DC: American Society for Microbiology; 1986:226-235.

10. Goyert SM, Ferrero E. Biochemical analysis of myeloid antigen and cDNA expression of gp55 (CD14). In: McMichael AJ, ed. *Leucocyte Typing III: White Cell Differentiation Antigens*. New York, NY: Oxford University Press; 1987:613-619.
11. van de Winkel JGJ, Capel JA. Human IgG Fc receptor heterogeneity: molecular aspects and clinical implication. *Immunol Today*. 1993;14:215-221.
12. van de Winkel JGJ, Anderson CL. Biology of human immunoglobulin G Fc receptors. *J Leukoc Biol*. 1991;49:511-524.
13. Bernstein ID, Self S. Joint report of the myeloid section of the Second International Workshop on Human Leukocyte Differentiation Antigens. In: Reinherz EL, Haynes BF, Nadler LM, Bernstein ID, eds. *Leukocyte Typing II: Human Myeloid and Hematopoietic Cells*. New York, NY: Springer-Verlag; 1986;3:1-25.
14. Jayaram Y, Hogg N. Surface expression of CD14 molecules on human neutrophils. In: Knapp W, Dörken B, Gilks W, et al, eds. *Leucocyte Typing IV: White Cell Differentiation Antigens*. New York, NY: Oxford University Press; 1989:796-797.
15. Shen L, Guyre PM, Fanger MW. Polymorphonuclear leukocyte function triggered through the high affinity Fc receptor for monomeric IgG. *J Immunol*. 1987;139:534-548.
16. Petroni KC, Shen L, Guyre PM. Modulation of human polymorphonuclear leukocyte IgG Fc receptors and Fc receptor-mediated functions by IFN- $\gamma$  and glucocorticoids. *J Immunol*. 1988;140:3467-3472.
17. Buckle A-M, Hogg N. The effect of IFN- $\gamma$  and colony-stimulating factors on the expression of neutrophil cell membrane receptors. *J Immunol*. 1989;143:2295-2301.
18. Perussia B, Dayton ET, Lazarus R, Fanning V, Trinchieri G. Immune interferon induces the receptor for monomeric IgG<sub>1</sub> on human monocytic and myeloid cells. *J Exp Med*. 1983;158:1092-1113.
19. Fanger NA, Wardwell K, Shen L, Tedder TF, Guyre PM. Type I (CD64) and Type II (CD32) Fc $\gamma$  receptor-mediated phagocytosis by human blood dendritic cells. *J Immunol*. 1996;157:541-548.
20. Valerius T, Repp R, de Wit TPM, et al. Involvement of the high-affinity receptor for IgG (Fc $\gamma$ RI; CD64) in enhanced tumor cell cytotoxicity of neutrophils during granulocyte colony-stimulating factor therapy. *Blood*. 1993;82:931-939.
21. Repp R, Valerius T, Sandler A, et al. Neutrophils express the high affinity receptor for IgG (Fc $\gamma$ RI, CD64) after in vivo application of recombinant human granulocyte colony-stimulating factor. *Blood*. 1991;78:885-889.
22. Kerst JM, van de Winkel JGJ, Evans AH, et al. Granulocyte colony-stimulating factor induces hFc gamma RI (CD64 antigen)-positive neutrophils via an effect on myeloid precursor cells. *Blood*. 1993;6:1457-1464.
23. Ernst LK, van de Winkel JGJ, Chiu IM, Anderson CL. Three genes for the human high affinity Fc receptor for IgG (Fc gamma RI) encode four distinct transcription products. *J Biol Chem*. 1992;267:15692-15700.

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BD, логотип BD и другие товарные знаки являются собственностью компании Becton, Dickinson and Company. © 2010 BD