BD Oncomark CD15 FITC/CD117 PE

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7/2010

23-13734-01

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- Schwartz S, Heinecke A, Zimmermann M, et al. Expression of the C-kit receptor (CD117) is a feature of almost all subtypes of de novo acute myeloblastic leukemia (AML), including cytogenetically good-risk AML, and lacks prognostic significance. *Leuk Lymphoma*. 1999;34(1-2):85-94.
- Cascavilla N, Musto P, D'Arena G, et al. CD117 (c-kit) is a restricted antigen of acute myeloid leukemia and characterizes early differentiative levels of M5 FAB subtype. *Haematologica*. 1998;83(5):392-397.
- Macedo A, Orfao A, Gonzalez M, et al. Immunological detection of blast cell subpopulations in acute myeloblastic leukemia at diagnosis: implications for minimal residual disease studies. *Leukemia*. 1995;9:993-998.
- Wells SJ, Bray RA, Stempora LL, Farhi DC. CD117/CD34 expression in leukemic blasts. *Am J Clin Pathol*. 1996;106(2):192-195.
- Terstappen LWMM, Safford M, Könemann S, et al. Flow cytometric characterization of acute myeloid leukemia. Part II. Phenotypic heterogeneity at diagnosis. *Leukemia*. 1992;6:70-80.
- Cohen PL, Hoyer JD, Kurtin PJ, Dewald GW, Hanson CA. Acute myeloid leukemia with minimal differentiation. A multiple parameter study [see comments]. *Am J Clin Pathol*. 1998;109(1):32-38.
- Sun T, Sangaline R, Ryder J, et al. Gating strategy for immunophenotyping of leukemia and lymphoma. *Am J Clin Pathol*. 1997;108(2):152-157.
- Matsuo Y, Adachi T, Tsubota T, Imanishi J, Minowada J. Establishment and characterization of a novel megakaryoblastoid cell line, MOLM-1, from a patient with chronic myelogenous leukemia. *Human Cell*. 1991;4:261-264.
- Ashman LK, Cambareri A, Nguyen L, Bühring H-J. CD117 workshop panel report. In: Kishimoto T, Kikutani H, von dem Borne AEGK, et al, eds. *Leucocyte Typing VI: White Cell Differentiation Antigens*. New York, NY: Garland Publishing, Inc; 1997:816-818.
- Rothe G, Schmitz G. Consensus protocol for the flow cytometric immunophenotyping of hematopoietic malignancies. *Leukemia*. 1996;10:877-895.
- Stelzer GT, Marti G, Hurley A, McCoy P, Jr., Lovett EJ, Schwartz A. US-Canadian consensus recommendations on the immunophenotypic analysis of hematologic neoplasia by flow cytometry: standardization and validation of laboratory procedures. *Cytometry*. 1997;30:214-230.
- Protection of Laboratory Workers from Occupationally Acquired Infections—Second Edition; Approved Guideline*. Villanova, PA: National Committee for Clinical Laboratory Standards; 2001. NCCLS document M29-A2.
- Centers for Disease Control. Update: universal precautions for prevention of transmission of human immunodeficiency virus, hepatitis B virus, and other bloodborne pathogens in health-care settings. *MMWR*. 1988;37:377-388.
- Jackson AL, Warner NL. Preparation, staining, and analysis by flow cytometry of peripheral blood leukocytes. In: Rose NR, Friedman H, Fahey JL, eds. *Manual of Clinical Laboratory Immunology*. 3rd ed. Washington, DC: American Society for Microbiology; 1986:226-235.
- Skubitz K, Balke J, Ball E, et al. Report on the CD15 cluster workshop. In: Knapp W, Dörken B, Gilks W, et al, eds. *Leucocyte Typing IV: White Cell Differentiation Antigens*. New York, NY: Oxford University Press; 1989:800-805.
- Berardi AC, Wang A, Levine JD, Lopez P, Scadden DT. Functional isolation and characterization of human hematopoietic stem cells. *Science*. 1995;267:104-108.
- Gunji Y, Nakamura M, Osawa H, et al. Human primitive hematopoietic progenitor cells are more enriched in KIT^{low} cells than in KIT^{high} cells. *Blood*. 1993;82:3283-3289.
- Broudy VC, Lin N, Zsebo KM, et al. Isolation and characterization of a monoclonal antibody that recognizes the human c-kit receptor. *Blood*. 1992;79:338-346.
- Briddell RA, Broudy VC, Bruno E, Brandt JE, Srou EF, Hoffman R. Further phenotypic characterization and isolation of human hematopoietic progenitor cells using a monoclonal antibody to the c-kit receptor. *Blood*. 1992;79:3159-3167.
- Simmons P, Aylett G, Niutta S, To L, Juttner C, Ashman L. *c-kit* is expressed by primitive human hematopoietic cells that give rise to colony-forming cells in stroma-dependent or cytokine-supplemented culture. *Exp Hematol*. 1994;22:157-165.
- Olweus J, Terstappen LWMM, Thompson PA, Lund-Johansen F. Expression and function of receptors for stem cell factor and erythropoietin during lineage commitment of human hematopoietic progenitor cells. *Blood*. 1996;88:1594-1607.
- Yarden Y, Kuang W-J, Yang-Feng T, et al. Human proto-oncogene c-kit: a new cell surface receptor tyrosine kinase for an unidentified ligand. *EMBO J*. 1987;6:3341-3351.
- Hanjan SN, Kearney JF, Cooper MD. A monoclonal antibody (MMA) that identifies a differentiation antigen on human myelomonocytic cells. *Clin Immunol Immunopathol*. 1982;23:172-188.

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BD, логотип BD и другие товарные знаки являются собственностью компании Becton, Dickinson and Company. © 2010 BD