# BD Oncomark CD15 FITC/CD34 PE

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7/2010

23-13754-01

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1. Terstappen LWM, Safford M, Könemann S, et al. Flow cytometric characterization of acute myeloid leukemia. Part II. Phenotypic heterogeneity at diagnosis. *Leukemia*. 1992;6:70-80.
2. Reading CL, Estey EH, Huh YO, et al. Expression of unusual immunophenotype combinations in acute myelogenous leukemia. *Blood*. 1993;81(11):3083-3090.
3. Venditti A, Del Poeta G, Stasi R, et al. Triple immunofluorescence evaluation of CD15, CD34 and class II expression by flow cytometry in normal and leukemic bone marrows. *Haematologica*. 1993;78:359-363.
4. Macedo A, Orfao A, Gonzalez M, et al. Immunological detection of blast cell subpopulations in acute myeloblastic leukemia at diagnosis: implications for minimal residual disease studies. *Leukemia*. 1995;9:993-998.
5. Orfao A, Chillon MC, Bortoluci AM, et al. The flow cytometric pattern of CD34, CD15 and CD13 expression in acute myeloblastic leukemia is highly characteristic of the presence of PML-RARalpha gene rearrangements. *Haematologica*. 1999;84(5):405-412.
6. Lansdorp PM, Dougherty GJ, Humphries RK. CD34 epitopes. In: Knapp W, Dörken B, Gilks WR, et al, eds. *Leucocyte Typing IV: White Cell Differentiation Antigens*. New York, NY: Oxford University Press; 1989:826-827.
7. Rothe G, Schmitz G. Consensus protocol for the flow cytometric immunophenotyping of hematopoietic malignancies. *Leukemia*. 1996;10:877-895.
8. Stelzer GT, Marti G, Hurley A, McCoy P, Jr, Lovett EJ, Schwartz A. US-Canadian consensus recommendations on the immunophenotypic analysis of hematologic neoplasia by flow cytometry: standardization and validation of laboratory procedures. *Cytometry*. 1997;30:214-230.
9. *Protection of Laboratory Workers from Occupationally Acquired Infections—Second Edition; Approved Guideline*. Villanova, PA: National Committee for Clinical Laboratory Standards; 2001. NCCLS document M29-A2.
10. Centers for Disease Control. Update: universal precautions for prevention of transmission of human immunodeficiency virus, hepatitis B virus, and other bloodborne pathogens in health-care settings. *MMWR*. 1988;37:377-388.
11. Jackson AL, Warner NL. Preparation, staining, and analysis by flow cytometry of peripheral blood leukocytes. In: Rose NR, Friedman H, Fahey JL, eds. *Manual of Clinical Laboratory Immunology*. 3rd ed. Washington, DC: American Society for Microbiology; 1986:226-235.
12. Skubitz K, Balke J, Ball E, et al. Report on the CD15 cluster workshop. In: Knapp W, Dörken B, Gilks W, et al, eds. *Leucocyte Typing IV: White Cell Differentiation Antigens*. New York, NY: Oxford University Press; 1989:800-805.
13. Loken MR, Shah VO, Dattilio KL, Civin CI. Flow cytometric analysis of human bone marrow. II. Normal B-lymphocyte development. *Blood*. 1987;70:1316-1324.
14. Civin CI, Trischmann TM, Fackler MJ, et al. Report on the CD34 cluster workshop. In: Knapp W, Dörken B, Gilks WR, et al, eds. *Leucocyte Typing IV: White Cell Differentiation Antigens*. New York, NY: Oxford University Press; 1989:818-825.
15. Peschel C, Köller U. Cluster report: CD34. In: Knapp W, Dörken B, Gilks W, et al, eds. *Leucocyte Typing IV: White Cell Differentiation Antigens*. New York, NY: Oxford University Press; 1989:817-818.
16. Lansdorp P, Sutherland H, Eaves C. Selective expression of CD45 isoforms on functional subpopulations of CD34<sup>+</sup> hemopoietic cells from human bone marrow. *J Exp Med*. 1990;172:363-366.
17. Civin CI, Strauss LC, Brovall C, Fackler MJ, Schwartz JF, Shaper JH. Antigenic analysis of hematopoiesis. III. A hemopoietic progenitor cell surface antigen defined by a monoclonal antibody raised against KG-1a cells. *J Immunol*. 1984;133:157-165.
18. Hanjan SN, Kearney JF, Cooper MD. A monoclonal antibody (MMA) that identifies a differentiation antigen on human myelomonocytic cells. *Clin Immunol Immunopathol*. 1982;23:172-188.
19. Kurtzberg J, Denning SM, Nycum LM, Singer KH, Haynes BF. Immature human thymocytes can be driven to differentiate into nonlymphoid lineages by cytokines from thymic epithelial cells. *Proc Natl Acad Sci USA*. 1989;86:7575-7579.
20. Terstappen LWM, Huang S, Safford M, Lansdorp PM, Loken MR. Sequential generations of hematopoietic colonies derived from single nonlineage-committed CD34<sup>+</sup>CD38<sup>-</sup> progenitor cells. *Blood*. 1991;77:1218-1227.
21. Civin C, Banquerigo M, Strauss L, Loken MR. Antigenic analysis of hematopoiesis. VI. Characterization of MY-10-positive progenitor cells in normal human bone marrow. *Exp Hematol*. 1987;15:10-17.

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BD, логотип BD и другие товарные знаки являются собственностью компании Becton, Dickinson and Company. © 2010 BD