BD Oncomark CD5 FITC/CD10 PE/ CD19 PerCP-Cy5.5

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6/2010

23-13713-01

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- Loken MR, Shah VO, Dattilio KL, Civin CI. Flow cytometric analysis of human bone marrow. II. Normal B-lymphocyte development. *Blood*. 1987;70:1316-1324.
- Hoffkes HG, Schmidtke G, Uppenkamp M, Schmucker U. Multiparametric immunophenotyping of B cells in peripheral blood of healthy adults by flow cytometry. *Clin Diagn Lab Immunol*. 1996;3:30-36.
- Freedman AS, Nadler LM. Immunologic markers in non-Hodgkin's lymphoma. *Hematol Oncol Clin North Am*. 1991;5:871-889.
- Jennings CD, Foon KA. Recent advances in flow cytometry: application to the diagnosis of hematologic malignancy. *Blood*. 1997;90:2863-2892.
- Zola H. CD10 Workshop Panel report. In: Schlossman SF, Boumsell L, Gilks W, et al, eds. *Leucocyte Typing V: White Cell Differentiation Antigens*. New York, NY: Oxford University Press; 1995:505-507.
- Gadol N, Ault KA. Phenotypic and functional characterization of human Leu-1 (CD5) B cells. *Immunol Rev*. 1986;93:23.
- Schmidt CJ, Domenico L, Ward P, Barcos MP, Stewart CC. Aberrant antigen expression detected by multiparameter three color flow cytometry in intermediate and high grade B-cell lymphomas. *Leuk Lymphoma*. 1999;34:539-544.
- Cabezudo E, Carrara P, Morilla R, Matutes E. Quantitative analysis of CD79b, CD5 and CD19 in mature B-cell lymphoproliferative disorders. *Haematologica*. 1999;84:413-418.
- Almasri NM, Iturraspe JA, Braylan RC. CD10 expression in follicular lymphoma and large cell lymphoma is different from that of reactive lymph node follicles. *Arch Pathol Lab Med*. 1998;122:539-544.
- Pui CH, Rivera GK, Hancock ML, et al. Clinical significance of CD10 expression in childhood acute lymphoblastic leukemia. *Leukemia*. 1993;7:35-40.
- De Rossi G, Grossi C, Foa R, et al. Immunophenotype of acute lymphoblastic leukemia cells: the experience of the Italian Cooperative Group (Gimema). *Leuk Lymphoma*. 1993;9:221-228.
- Geisler CH, Larsen JK, E HN, et al. Prognostic importance of flow cytometric immunophenotyping of 540 consecutive patients with B-cell chronic lymphocytic leukemia. *Blood*. 1991;78:1795-1802.
- DiGiuseppe JA, Borowitz MJ. Clinical utility of flow cytometry in the chronic lymphoid leukemias. *Semin Oncol*. 1998;25:6-10.

14. Baldini L, Cro L, Cortelezzi A, et al. Immunophenotypes in "classical" B-cell chronic lymphocytic leukemia. Correlation with normal cellular counterpart and clinical findings. *Cancer*. 1990;66(8):1738-1742.
15. Spier CM, Grogan TM, Fielder K, Richter L, Rangel C. Immunophenotypes in "well-differentiated" lymphoproliferative disorders, with emphasis on small lymphocytic lymphoma. *Hum Pathol*. 1986;17(11):1126-1136.
16. Engleman EG, Warnke R, Fox RI, Dille J, Benike CJ, Levy R. Studies of a human T lymphocyte antigen recognized by a monoclonal antibody. *Proc Natl Acad Sci USA*. 1981;78:1791-1795.
17. Nadler LM. B Cell/Leukemia Panel Workshop: summary and comments. In: Reinherz EL, Haynes BF, Nadler LM, Bernstein ID, eds. *Leukocyte Typing II: Human B Lymphocytes*. New York, NY: Springer-Verlag; 1986;2:3-43.
18. Rothe G, Schmitz G. Consensus protocol for the flow cytometric immunophenotyping of hematopoietic malignancies. *Leukemia*. 1996;10:877-895.
19. Stelzer GT, Marti G, Hurley A, McCoy P, Jr., Lovett EJ, Schwartz A. US-Canadian consensus recommendations on the immunophenotypic analysis of hematologic neoplasia by flow cytometry: standardization and validation of laboratory procedures. *Cytometry*. 1997;30:214-230.
20. *Protection of Laboratory Workers from Occupationally Acquired Infections—Second Edition; Approved Guideline*. Villanova, PA: National Committee for Clinical Laboratory Standards; 2001. NCCLS document M29-A2.
21. Centers for Disease Control. Update: universal precautions for prevention of transmission of human immunodeficiency virus, hepatitis B virus, and other bloodborne pathogens in health-care settings. *MMWR*. 1988;37:377-388.
22. Jackson AL, Warner NL. Preparation, staining, and analysis by flow cytometry of peripheral blood leukocytes. In: Rose NR, Friedman H, Fahey JL, eds. *Manual of Clinical Laboratory Immunology*. 3rd ed. Washington, DC: American Society for Microbiology; 1986:226-235.
23. Knowles RW. Immunocytochemical analysis of the T-cell-specific antigens. In: Reinherz EL, Haynes BF, Nadler LM, Bernstein ID, eds. *Leukocyte Typing II: Human T Lymphocytes*. New York, NY: Springer-Verlag; 1986;1:259-288.
24. Greaves MF. Monoclonal antibodies as probes for leukemic heterogeneity and hematopoietic differentiation. In: Knapp W, ed. *Leukemia Markers*. New York, NY: Academic Press; 1981:19.
25. Letarte M, Vera S, Tran R, et al. Common acute lymphocytic leukemia antigen is identical to neutral endopeptidase. *J Exp Med*. 1988;168:1247-1253.
26. Moldenhauer G, Dörken B, Schwartz R, Pezzutto A, Knops J, Hammerling GJ. Analysis of ten B lymphocyte-specific workshop monoclonal antibodies. In: Reinherz EL, Haynes BF, Nadler LM, Bernstein ID, eds. *Leukocyte Typing II: Human B Lymphocytes*. New York, NY: Springer-Verlag; 1986;2:61-67.
27. Ledbetter JA, Evans RL, Lipinski M, Cunningham-Rundles C, Good RA, Herzenberg LA. Evolutionary conservation of surface molecules that distinguish T-lymphocyte helper/inducer and T cytotoxic/suppressor subpopulations in mouse and man. *J Exp Med*. 1981;153:310-323.
28. Ledbetter JA, Frankel AE, Herzenberg LA. Human Leu T-cell differentiation antigens: quantitative expression on normal lymphoid cells and cell lines. In: Hammerling G, Hammerling U, Kearney J, eds. *Monoclonal Antibodies and T Cell Hybridomas: Perspectives and Technical Notes*. New York, NY: Elsevier/North Holland; 1981:16-22.
29. Warnke RA, Levy R. Detection of T and B antigens with hybridoma monoclonal antibodies: a biotin-avidin-horseradish peroxidase method. *J Histochem Cytochem*. 1980;28:771-776.
30. Warnke R, Miller R, Grogan T, Pederson M, Dille J, Levy R. Immunologic phenotype in 30 patients with diffuse large-cell lymphoma. *N Eng J Med*. 1980;303:293-300.
31. Zipf RF, Fox R, Dille J, Levy R. Definition of the high risk ALL patient by immunologic phenotyping with monoclonal antibodies. *Cancer Res*. 1981;41:4786.
32. Royston I, Majda JA, Baird SM, Meserve BL, Griffiths JC. Human T-cell antigens defined by monoclonal antibodies: the 65,000-dalton antigen of T cells (T65) is also found on chronic lymphocytic leukemia cells bearing surface immunoglobulin. *J Immunol*. 1980;125:725.
33. Le Bien TW, McCormack RT. The common acute lymphoblastic leukemia antigen (CD10): emancipation from a functional enigma. *Blood*. 1989;73:625-635.
34. Reichert T, DeBruyère M, Deneys V, et al. Lymphocyte subset reference ranges in adult Caucasians. *Clin Immunol Immunopath*. 1991;60:190-208.
35. Tedder T, Zhou L-J, Engel P. The CD19/CD21 signal transduction complex of B lymphocytes. *Immunol Today*. 1994;15:437-442.

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BD, логотип BD и другие товарные знаки являются собственностью компании Becton, Dickinson and Company. © 2010 BD