BD Oncomark CD7 FITC/CD33 PE

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7/2010

23-13733-01

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1. Terstappen LWM, Safford M, Könemann S, et al. Flow cytometric characterization of acute myeloid leukemia. Part II. Phenotypic heterogeneity at diagnosis. *Leukemia*. 1992;6:70-80.
2. Del Poeta G, Stasi R, Venditti A, et al. CD7 expression in acute myeloid leukemia. *Leuk Lymphoma*. 1995; 17(1-2):111-119.
3. Macedo A, Orfao A, Vidriales MB, et al. Characterization of aberrant phenotypes in acute myeloblastic leukemia. *Ann Hematol*. 1995;70(4):189-194.
4. Kornblau SM, Thall P, Huh YO, Estey E, Andreeff M. Analysis of CD7 expression in acute myelogenous leukemia: martingale residual plots combined with "optimal" cutpoint analysis reveals absence of prognostic significance. *Leukemia*. 1995;9(10):1735-1741.
5. Tien HF, Chou CC, Wang CH, Chang CH, Hsing CC. Putative normal counterparts of leukaemic cells from CD7-positive acute myeloid leukaemia can be demonstrated in human haemopoietic tissues. *Br J Haematol*. 1996;94(3):501-506.
6. Saxena A, Sheridan DP, Card RT, McPeck AM, Mewdell CC, Skinnider LF. Biologic and clinical significance of CD7 expression in acute myeloid leukemia. *Am J Hematol*. 1998;58(4):278-284.
7. Macedo A, Orfao A, Vidriales MB, et al. Characterization of aberrant phenotypes in acute myeloblastic leukemia. *Ann Hematol*. 1995;70:189-194.
8. Reiter C. Cluster report: CD7. In: Knapp W, Dörken B, Gils WR, et al, eds. *Leucocyte Typing IV: White Cell Differentiation Antigens*. New York, NY: Oxford University Press; 1989:341-342.
9. Rothe G, Schmitz G. Consensus protocol for the flow cytometric immunophenotyping of hematopoietic malignancies. *Leukemia*. 1996;10:877-895.
10. Stelzer GT, Marti G, Hurley A, McCoy P, Jr., Lovett EJ, Schwartz A. US-Canadian consensus recommendations on the immunophenotypic analysis of hematologic neoplasia by flow cytometry: standardization and validation of laboratory procedures. *Cytometry*. 1997;30:214-230.
11. *Protection of Laboratory Workers from Occupationally Acquired Infections—Second Edition; Approved Guideline*. Villanova, PA: National Committee for Clinical Laboratory Standards; 2001. NCCLS document M29-A2.
12. Centers for Disease Control. Update: universal precautions for prevention of transmission of human immunodeficiency virus, hepatitis B virus, and other bloodborne pathogens in health-care settings. *MMWR*. 1988;37:377-388.

13. Jackson AL, Warner NL. Preparation, staining, and analysis by flow cytometry of peripheral blood leukocytes. In: Rose NR, Friedman H, Fahey JL, eds. *Manual of Clinical Laboratory Immunology*. 3rd ed. Washington, DC: American Society for Microbiology; 1986:226-235.
14. Link M, Warnke R, Finlay J, et al. A single monoclonal antibody identifies T-cell lineage of childhood lymphoid malignancies. *Blood*. 1983;62:722-728.
15. Palker TJ, Searce RM, Hensley LL, Ho W, Haynes BF. Comparison of the CD7 (3A1) group of T cell workshop antibodies. In: Reinherz EL, Haynes BF, Nadler LM, Bernstein ID, eds. *Leukocyte Typing II: Human T Lymphocytes*. New York, NY: Springer-Verlag; 1986;1:303-313.
16. Köller U, Peschel CH. Cluster report: CD33. In: Knapp W, Dörken B, Gilks WR, et al, eds. *Leucocyte Typing IV: White Cell Differentiation Antigens*. New York, NY: Oxford University Press; 1989:812-813.
17. Weiss LM, Crabtree GS, Rouse RV, Warnke RA. Morphologic and immunologic characterization of 50 peripheral T-cell lymphomas. *Am J Pathol*. 1985;118:316-324.
18. Terstappen LWMM, Hollander Z, Meiners H, Loken MR. Quantitative comparison of myeloid antigens on five lineages of mature peripheral blood cells. *J Leuk Biol*. 1990;48:138-148.
19. Foon KA, Todd III RF. Immunologic classification of leukemia and lymphoma. *Blood*. 1986;68:1-31.
20. Bernstein ID, Singer JW, Andrews RG, et al. Treatment of acute myeloid leukemia cells *in vitro* with a monoclonal antibody recognizing a myeloid differentiation antigen allows normal progenitor cells to be expressed. *J Clin Invest*. 1987;79:1153-1159.

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BD, логотип BD и другие товарные знаки являются собственностью компании Becton, Dickinson and Company. © 2010 BD