# BD Oncomark FMC7 FITC/CD23 PE/CD19 PerCP-Cy5.5

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6/2010

23-13714-01

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1. Tworek JA, Singleton TP, Schnitzer B, Hsi ED, Ross CW. Flow cytometric and immunohistochemical analysis of small lymphocytic lymphoma, mantle cell lymphoma, and plasmacytoid small lymphocytic lymphoma. *Am J Clin Pathol*. 1998;110:582-589.
2. Hoffkes HG, Schmidke G, Uppenkamp M, Schmucker U. Multiparametric immunophenotyping of B cells in peripheral blood of healthy adults by flow cytometry. *Clin Diagn Lab Immunol*. 1996;3:30-36.
3. Huh YO, Pugh WC, Kanjarjian HM, et al. Detection of subgroups of chronic B-cell leukemias by FMC7 monoclonal antibody. *Am J Clin Pathol*. 1994;101:283-289.
4. DiGiuseppe JA, Borowitz MJ. Clinical utility of flow cytometry in the chronic lymphoid leukemias. *Semin Oncol*. 1998;25:6-10.
5. Hübl W, Iturraspe J, Braylan RC. FMC7 antigen expression on normal and malignant B-cells can be predicted by expression of CD20. *Cytometry*. 1998;34:71-74.
6. Kilo MN, Dorfman DM. The utility of flow cytometric immunophenotypic analysis in the distinction of small lymphocytic lymphoma/chronic lymphocytic leukemia from mantle cell lymphoma. *Am J Clin Pathol*. 1996;105:451-457.
7. Geisler CH, Larsen JK, E HN, et al. Prognostic importance of flow cytometric immunophenotyping of 540 consecutive patients with B-cell chronic lymphocytic leukemia. *Blood*. 1991;78:1795-1802.
8. Marti GE, Zenger V, Caproaso NE, et al. Antigenic expression of B-cell chronic lymphocytic leukemic lymphocytes. *Anal Quant Cytol Histol*. 1989;11:315-323.
9. Matutes E, Owusu-Ankomah K, Morilla R, et al. The immunological profile of B-cell disorders and proposal of a scoring system for the diagnosis of CLL. *Leukemia*. 1994;8(10):1640-1645.
10. Brooks DA, Beckman IGR, Bradley J, McNamara PJ, Thomas ME, Zola H. Human lymphocyte markers defined by antibodies derived from somatic cell hybrids. IV. A monoclonal antibody reacting specifically with a subpopulation of human B lymphocytes. *J Immunol*. 1981;126:1373-1377.
11. Kintner C, Sugden B. Identification of antigenic determinants unique to the surfaces of cells transformed by Epstein-Barr Virus. *Nature*. 1981;294:458-460.
12. Nadler LM. B Cell/Leukemia Panel Workshop: summary and comments. In: Reinherz EL, Haynes BF, Nadler LM, Bernstein ID, eds. *Leukocyte Typing II: Human B Lymphocytes*. New York, NY: Springer-Verlag; 1986;2:3-43.
13. Rothe G, Schmitz G. Consensus protocol for the flow cytometric immunophenotyping of hematopoietic malignancies. *Leukemia*. 1996;10:877-895.
14. Stelzer GT, Marti G, Hurley A, McCoy P, Jr., Lovett EJ, Schwartz A. US-Canadian consensus recommendations on the immunophenotypic analysis of hematologic neoplasia by flow cytometry: standardization and validation of laboratory procedures. *Cytometry*. 1997;30:214-230.
15. *Protection of Laboratory Workers from Occupationally Acquired Infections—Second Edition; Approved Guideline*. Villanova, PA: National Committee for Clinical Laboratory Standards; 2001. NCCLS document M29-A2.
16. Centers for Disease Control. Update: universal precautions for prevention of transmission of human immunodeficiency virus, hepatitis B virus, and other bloodborne pathogens in health-care settings. *MMWR*. 1988;37:377-388.
17. Jackson AL, Warner NL. Preparation, staining, and analysis by flow cytometry of peripheral blood leukocytes. In: Rose NR, Friedman H, Fahey JL, eds. *Manual of Clinical Laboratory Immunology*. 3rd ed. Washington, DC: American Society for Microbiology; 1986:226-235.
18. Nadler LM. Cluster report: CD23. In: Reinherz EL, Haynes BF, Nadler LM, Bernstein ID, eds. *Leukocyte Typing II: Human B Lymphocytes*. New York, NY: Springer-Verlag; 1986;2:25-26.
19. Gordon J, Webb A, Walker L, Guy G, Rowe M. Evidence of an association between CD23 and the receptor for a low molecular weight B-cell growth factor. *Eur J Immunol*. 1986;16:1627.
20. Yukawa K, Kikutani H, Owaki H, et al. A B-cell-specific differentiation antigen, CD23, is a receptor for IgE (Fc  $\epsilon$  R) on lymphocytes. *J Immunol*. 1987;138:2576-2580.
21. Thorley-Lawson D, Nadler L, Bhan A, Schooley R. Blast-2 (EBVCS), an early cell surface marker of human B cell activation, is superinduced by Epstein Barr Virus. *J Immunol*. 1985;134:3007-3012.
22. Moldenhauer G, Dörken B, Schwartz R, Pezzutto A, Knops J, Hammerling GJ. Analysis of ten B lymphocyte-specific workshop monoclonal antibodies. In: Reinherz EL, Haynes BF, Nadler LM, Bernstein ID, eds. *Leukocyte Typing II: Human B Lymphocytes*. New York, NY: Springer-Verlag; 1986;2:61-67.
23. Bloem AC, Chand MA, Dollekamp I, Rijkers GT. Functional properties of human B cell subpopulations defined by monoclonal antibodies HB4 and FMC7. *J Immunol*. 1988;140:768-773.
24. Hassan IB, Hagberg H, Sundström C. Immunophenotype of hairy-cell leukemia. *Eur J Haematol*. 1990;45:172-176.

25. Zola H, Neoh SH, Potter A, Melo JV, De Oliveria MSP, Catovsky D. Markers of differentiated B cell leukaemia: CD22 antibodies and FMC7 react with different molecules. *Disease Markers*. 1987;5:227-235.
26. Matutes E, Morilla R, Owusu-Ankomah K, Houlihan A, Meeus P, Catovsky D. The immunophenotype of hairy cell leukemia (HCL). Proposal for a scoring system to distinguish HCL from B-cell disorders with hairy or villous lymphocytes. *Leuk Lymphoma*. 1994;14:57-61.
27. Kikutani H, Inui S, Sato R, et al. Molecular structure of human lymphocyte receptor for immunoglobulin E. *Cell*. 1986;47:657-665.
28. Gordon J, Rowe M, Walker L, Guy G. Ligation of the CD23,p45 (Blast-2, EBVCS) antigen triggers the cell-cycle progression of activated B lymphocytes. *Eur J Immunol*. 1986;16:1075-1080.
29. Capron M, Jouault T, Prin L, et al. Functional study of a monoclonal antibody to IgE Fc receptor (FcεR2) of eosinophils, platelets, and macrophages. *J Exp Med*. 1986;164:72-89.
30. Reichert T, DeBruyère M, Deneys V, et al. Lymphocyte subset reference ranges in adult Caucasians. *Clin Immunol Immunopath*. 1991;60:190-208.
31. Tedder T, Zhou L-J, Engel P. The CD19/CD21 signal transduction complex of B lymphocytes. *Immunol Today*. 1994;15:437-442.
32. Loken MR, Shah VO, Dattilio KL, Civin CI. Flow cytometric analysis of human bone marrow. II. Normal B-lymphocyte development. *Blood*. 1987;70:1316-1324.

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BD, логотип BD и другие товарные знаки являются собственностью компании Becton, Dickinson and Company. © 2010 BD