BD Oncomark Anti-HLA-DR FITC/CD13 PE

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7/2010

23-13736-01

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- Guglielmi C, Martelli MP, Diverio D, et al. Immunophenotype of adult and childhood acute promyelocytic leukaemia: correlation with morphology, type of PML gene breakpoint and clinical outcome. A cooperative Italian study on 196 cases. *Br J Haematol*. 1998;102(4):1035-1041.
- Paietta E, Andersen J, Gallagher R, et al. The immunophenotype of acute promyelocytic leukemia (APL): an ECOG study. *Leukemia*. 1994;8:1108-1112.
- Macedo A, Orfao A, Vidriales MB, et al. Characterization of aberrant phenotypes in acute myeloblastic leukemia. *Ann Hematol*. 1995;70(4):189-194.
- Terstappen LWM, Safford M, Könemann S, et al. Flow cytometric characterization of acute myeloid leukemia. Part II. Phenotypic heterogeneity at diagnosis. *Leukemia*. 1992;6:70-80.
- Lampson LA, Levy R. Two populations of Ia-like molecules on a human B cell line. *J Immunol*. 1980;125:293-299.
- Rothe G, Schmitz G. Consensus protocol for the flow cytometric immunophenotyping of hematopoietic malignancies. *Leukemia*. 1996;10:877-895.
- Stelzer GT, Marti G, Hurley A, McCoy P, Jr, Lovett EJ, Schwartz A. US-Canadian consensus recommendations on the immunophenotypic analysis of hematologic neoplasia by flow cytometry: standardization and validation of laboratory procedures. *Cytometry*. 1997;30:214-230.
- Protection of Laboratory Workers from Occupationally Acquired Infections—Second Edition; Approved Guideline*. Villanova, PA: National Committee for Clinical Laboratory Standards; 2001. NCCLS document M29-A2.
- Centers for Disease Control. Update: universal precautions for prevention of transmission of human immunodeficiency virus, hepatitis B virus, and other bloodborne pathogens in health-care settings. *MMWR*. 1988;37:377-388.
- Jackson AL, Warner NL. Preparation, staining, and analysis by flow cytometry of peripheral blood leukocytes. In: Rose NR, Friedman H, Fahey JL, eds. *Manual of Clinical Laboratory Immunology*. 3rd ed. Washington, DC: American Society for Microbiology; 1986:226-235.
- Brodsky FM. A matrix approach to human class II histocompatibility antigens: reactions of four monoclonal antibodies with the products of nine haplotypes. *Immunogenetics*. 1984;19:179-194.
- Gadd S. Cluster report: CD13. In: Knapp W, Dörken B, Gilks WR, et al, eds. *Leucocyte Typing IV: White Cell Differentiation Antigens*. New York, NY: Oxford University Press; 1989:782-784.
- Look A, Ashmun R, Shapiro L, et al. Report on the CD13 (aminopeptidase N) cluster workshop. In: Knapp W, Dörken B, Gilks W, et al, eds. *Leucocyte Typing IV: White Cell Differentiation Antigens*. New York, NY: Oxford University Press; 1989:784-787.
- Tomkinson B, Wagner D, Nelson D, Sullivan J. Activated lymphocytes during acute Epstein-Barr virus infection. *J Immunol*. 1987;139:3802-3807.
- Levacher M, Tallet S, Dazza M, Dournon E, Rouveix B, Pocard J. T activation marker evaluation in ARC patients treated with AZT: comparison with CD4⁺ lymphocyte count in non-progressors and progressors towards AIDS. *Clin Exp Immunol*. 1990;81:177-182.
- Stites DP, Casavant CH, McHugh TM, et al. Flow cytometric analysis of lymphocyte phenotypes in AIDS using monoclonal antibodies and simultaneous dual immunofluorescence. *Clin Immunol Immunopathol*. 1986;38:161-177.
- Terstappen LWM, Hollander Z, Meiners H, Loken MR. Quantitative comparison of myeloid antigens on five lineages of mature peripheral blood cells. *J Leuk Biol*. 1990;48:138-148.
- Warnke RA, Levy R. Detection of T and B antigens with hybridoma monoclonal antibodies: a biotin-avidin-horseradish peroxidase method. *J Histochem Cytochem*. 1980;28:771-776.
- Edwards JA, Durant BM, Jones DB, Evans PR, Smith JL. Differential expression of HLA class II antigens in fetal human spleen: relationship of HLA-DP, DQ, and DR to immunoglobulin expression. *J Immunol*. 1985;137:490-497.
- Warnke R, Miller R, Grogan T, Pederson M, Dille J, Levy R. Immunologic phenotype in 30 patients with diffuse large-cell lymphoma. *N Engl J Med*. 1980;303:293-300.
- Engleman EG, Warnke R, Fox RI, Dille J, Benike CJ, Levy R. Studies of a human T lymphocyte antigen recognized by a monoclonal antibody. *Proc Natl Acad Sci USA*. 1981;78:1791-1795.
- Zipf RF, Fox R, Dille J, Levy R. Definition of the high risk ALL patient by immunologic phenotyping with monoclonal antibodies. *Cancer Res*. 1981;41:4786.
- Foon KA, Todd III RF. Immunologic classification of leukemia and lymphoma. *Blood*. 1986;68:1-31.
- Greaves MF, Chan LC, Furely AJW, Watt SM, Molgaard HV. Lineage promiscuity in hematopoietic differentiation and leukemia. *Blood*. 1986;67:1-11.
- Hurwitz CA, Loken MR, Graham ML, et al. Asynchronous antigen expression in B lineage acute lymphoblastic leukemia. *Blood*. 1988;72:299-307.
- Greaves MF, Harii G, Newman RA, Sutherland DR, Ritter MA, Ritz J. Selective expression of the common acute lymphoblastic leukemia (gp100) antigen on immature lymphoid cells and their malignant counterparts. *Blood*. 1983;61:628-639.

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BD, логотип BD и другие товарные знаки являются собственностью компании Becton, Dickinson and Company. © 2010 BD