BD Simultest CD5 FITC/CD19 PE

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7/2010

23-13775-02

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- Bain B, Barnett D, Linch D, Matutes E, Reilly J. Revised guideline on immunophenotyping in acute leukaemias and chronic lymphoproliferative disorders. *Clin Lab Haem.* 2002;24:1-13.
- Communications in Clinical Cytometry. *Cytometry.* 2001;46:23-27.
- Engleman EG, Warnke R, Fox RI, Dilley J, Benike CJ, Levy R. Studies of a human T lymphocyte antigen recognized by a monoclonal antibody. *Proc Natl Acad Sci USA.* 1981;78:1791-1795.
- Nadler LM. B Cell/Leukemia Panel Workshop: summary and comments. In: Reinherz EL, Haynes BF, Nadler LM, Bernstein ID, eds. *Leukocyte Typing II: Human B Lymphocytes.* New York, NY: Springer-Verlag; 1986;2:3-43.
- Rothe G, Schmitz G. Consensus protocol for the flow cytometric immunophenotyping of hematopoietic malignancies. *Leukemia.* 1996;10:877-895.
- Stelzer GT, Marti G, Hurley A, McCoy P, Jr, Lovett EJ, Schwartz A. US-Canadian consensus recommendations on the immunophenotypic analysis of hematologic neoplasia by flow cytometry: standardization and validation of laboratory procedures. *Cytometry.* 1997;30:214-230.
- Protection of Laboratory Workers from Occupationally Acquired Infections—Second Edition; Approved Guideline.* Villanova, PA: National Committee for Clinical Laboratory Standards; 2001. NCCLS document M29-A2.
- Centers for Disease Control. Update: universal precautions for prevention of transmission of human immunodeficiency virus, hepatitis B virus, and other bloodborne pathogens in health-care settings. *MMWR.* 1988;37:377-388.
- Jackson AL, Warner NL. Preparation, staining, and analysis by flow cytometry of peripheral blood leukocytes. In: Rose NR, Friedman H, Fahey JL, eds. *Manual of Clinical Laboratory Immunology.* 3rd ed. Washington, DC: American Society for Microbiology; 1986:226-235.
- Knowles RW. Immunochemical analysis of the T-cell-specific antigens. In: Reinherz EL, Haynes BF, Nadler LM, Bernstein ID, eds. *Leukocyte Typing II: Human T Lymphocytes.* New York, NY: Springer-Verlag; 1986;1:259-288.
- Moldenhauer G, Dörken B, Schwartz R, Pezzutto A, Knops J, Hammerling GJ. Analysis of ten B lymphocyte-specific workshop monoclonal antibodies. In: Reinherz EL, Haynes BF, Nadler LM, Bernstein ID, eds. *Leukocyte Typing II: Human B Lymphocytes.* New York, NY: Springer-Verlag; 1986;2:61-67.
- Ledbetter JA, Evans RL, Lipinski M, Cunningham-Rundles C, Good RA, Herzenberg LA. Evolutionary conservation of surface molecules that distinguish T-lymphocyte helper/inducer and T cytotoxic/suppressor subpopulations in mouse and man. *J Exp Med.* 1981;153:310-323.
- Ledbetter JA, Frankel AE, Herzenberg LA. Human Leu T-cell differentiation antigens: quantitative expression on normal lymphoid cells and cell lines. In: Hämmerling G, Hämmerling U, Kearney J, eds. *Monoclonal Antibodies and T Cell Hybridomas: Perspectives and Technical Notes.* New York, NY: Elsevier/North Holland; 1981:16-22.
- Warnke RA, Levy R. Detection of T and B antigens with hybridoma monoclonal antibodies: a biotin-avidin-horseradish peroxidase method. *J Histochem Cytochem.* 1980;28:771-776.
- Warnke R, Miller R, Grogan T, Pederson M, Dilley J, Levy R. Immunologic phenotype in 30 patients with diffuse large-cell lymphoma. *N Eng J Med.* 1980;303:293-300.
- Zipf RF, Fox R, Dilley J, Levy R. Definition of the high risk ALL patient by immunologic phenotyping with monoclonal antibodies. *Cancer Res.* 1981;41:4786.
- Gadol N, Ault KA. Phenotypic and functional characterization of human Leu-1 (CD5) B cells. *Immunol Rev.* 1986;93:23.
- Casali P, Burastero SE, Nakamura M, Inghirami G, Notkins AL. Human lymphocytes making rheumatoid factor and antibody to ssDNA belong to Leu-1⁺ B-cell subset. *Science.* 1987;236:77-80.
- Royston I, Majda JA, Baird SM, Meserve BL, Griffiths JC. Human T-cell antigens defined by monoclonal antibodies: the 65,000-dalton antigen of T cells (T65) is also found on chronic lymphocytic leukemia cells bearing surface immunoglobulin. *J Immunol.* 1980;125:725.
- Le Bien TW, McCormack RT. The common acute lymphoblastic leukemia antigen (CD10): emancipation from a functional enigma. *Blood.* 1989;73:625-635.
- Reichert T, DeBruyère M, Deneyns V, et al. Lymphocyte subset reference ranges in adult Caucasians. *Clin Immunol Immunopath.* 1991;60:190-208.
- Tedder T, Zhou L-J, Engel P. The CD19/CD21 signal transduction complex of B lymphocytes. *Immunol Today.* 1994;15:437-442.
- Loken MR, Shah VO, Dattilio KL, Civin CI. Flow cytometric analysis of human bone marrow. II. Normal B-lymphocyte development. *Blood.* 1987;70:1316-1324.

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BD, логотип BD и другие товарные знаки являются собственностью компании Becton, Dickinson and Company. © 2010 BD