TriTest CD4 FITC/ CD8 PE/CD3 PerCP Reagent

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1. Schmidt RE. Monoclonal antibodies for diagnosis of immunodeficiencies. *Blut*. 1989;59:200-206.
2. Nicholson JKA. Use of flow cytometry in the evaluation and diagnosis of primary and secondary immunodeficiency diseases. *Arch Pathol Lab Med*. 1989;113:598-605.
3. Foucar K, Goeken JA. Clinical application of immunologic techniques to the diagnosis of lymphoproliferative and immunodeficiency disorders. *Lab Med*. 1982;13:403-413.
4. Cohen SB, Weetman AP. Activated interstitial and intraepithelial thyroid lymphocytes in autoimmune thyroid disease. *Acta Endocrinol*. 1988;119:161-166.

5. Smolen JS, Chused TM, Leiserson WM, Reeves JP, Alling D, Steinberg AD. Heterogeneity of immunoregulatory T-cell subsets in systemic lupus erythematosus: correlation with clinical features. *Am J Med*. 1982;72:783-790.
6. Giorgi JV, Hultin LE. Lymphocyte subset alterations and immunophenotyping by flow cytometry in HIV disease. *Clin Immunol Newslett*. 1990;10:55-61.
7. Landay A, Ohlsson-Wilhelm B, Giorgi JV. Application of flow cytometry to the study of HIV infection. *AIDS*. 1990;4:479-497.
8. Antel J, Bania M, Noronha A, Neely S. Defective suppressor cell function mediated by T8⁺ cell lines from patients with progressive multiple sclerosis. *J Immunol*. 1986;137:3436-3439.
9. Gratama JW, Naipal A, Olijans P, et al. T-lymphocyte repopulation and differentiation after bone marrow transplantation: early shifts in the ratio between T4⁺ and T8⁺ T lymphocytes correlate with the occurrence of acute graft-versus-host disease. *Blood*. 1984;63:1416-1423.
10. Bishop GA, Hall BM, Duggin GG, Horvath JS, Sheil AGR, Tiller DJ. Immunopathology of renal allograft rejection analyzed with monoclonal antibodies to mononuclear cell markers. *Kidney Internat*. 1986;29:708-717.
11. Wolde-Mariam W, Peter JB. Recent diagnostic advances in cellular immunology. *Diag Med*. 1984;7:25-32.
12. Centers for Disease Control. 1997 Revised guidelines for performing CD4⁺ T-cell determinations in persons with human immunodeficiency virus (HIV). *MMWR*. 1997;46(No. RR-2):1-29.
13. Nicholson JKA, Jones BM, Hubbard M. CD4 T-lymphocyte determinations on whole blood specimens using a single-tube three-color assay. *Cytometry*. 1993;14:685-689.
14. Nicholson J, Kidd P, Mandy F, Livnat D, Kagan J. Three-color supplement to the NIAID DAIDS guideline for flow cytometric immunophenotyping. *Cytometry*. 1996;26:227-230.
15. Nicholson JKA, Hubbard M, Jones BM. Use of CD45 fluorescence and side-scatter characteristics for gating lymphocytes when using the whole blood lysis procedure and flow cytometry. *Cytometry*. 1996;26:16-21.
16. Evans RL, Wall DW, Platsoucas CD, et al. Thymus-dependent membrane antigens in man: inhibition of cell-mediated lympholysis by monoclonal antibodies to the T_{H2} antigen. *Proc Natl Acad Sci USA*. 1981;78:544-548.

17. Bernard A, Boumsell L, Hill C. Joint report of the first international workshop on human leucocyte differentiation antigens by the investigators of the participating laboratories: T2 protocol. In: Bernard A, Boumsell L, Dausset J, Milstein C, Schlossman SF, eds. *Leukocyte Typing*. New York, NY: Springer-Verlag; 1984:25-60.
18. Wood GS, Warner NL, Warnke RA. Anti-Leu-3/T4 antibodies react with cells of monocyte/macrophage and Langerhans lineage. *J Immunol*. 1983;131:212-216.
19. Haynes BF. Summary of T-cell studies performed during the Second International Workshop and Conference on Human Leukocyte Differentiation Antigens. In: Reinherz EL, Haynes BF, Nadler LM, Bernstein ID, eds. *Leukocyte Typing II: Human T Lymphocytes*. New York, NY: Springer-Verlag; 1986:3-30.
20. Kan EAR, Wang CY, Wang LC, Evans RL. Noncovalently bonded subunits of 22 and 28 kd are rapidly internalized by T cells reacted with Anti-Leu-4 antibody. *J Immunol*. 1983;131:536-539.
21. Knowles RW. Immunochemical analysis of the T-cell-specific antigens. In: Reinherz EL, Haynes BF, Nadler LM, Bernstein ID, eds. *Leukocyte Typing II: Human T Lymphocytes*. New York, NY: Springer-Verlag; 1986;1:259-288.
22. Organizing Committee of the Fourth International Workshop on Human Leukocyte Differentiation Antigens. Appendix A: CD guide. In: Knapp W, Dörken B, Gilks WR, et al, eds. *Leukocyte Typing IV: White Cell Differentiation Antigens*. New York, NY: Oxford University Press Inc; 1989:1074-1093.
23. Dalgleish AG, Beverley PCL, Clapham PR, Crawford DH, Greaves MF, Weiss RA. The CD4 (T4) antigen is an essential component of the receptor for the AIDS retrovirus. *Nature*. 1984;312:763-767.
24. Maddon PJ, Dalgleish AG, McDougal JS, Clapham PR, Weiss RA, Axel R. The T4 gene encodes the AIDS virus receptor and is expressed in the immune system and the brain. *Cell*. 1986;47:333-348.
25. Rudd CE, Burgess KE, Barber EK, Schlossman SF. Monoclonal antibodies to the CD4 and CD8 antigens precipitate variable amounts of CD4/CD8-associated p56^{lck} activity. In: Knapp W, Dörken B, Gilks WR, et al, eds. *Leukocyte Typing IV: White Cell Differentiation Antigens*. New York, NY: Oxford University Press; 1989:326-327.
26. Moebius U. Cluster report: CD8. In: Knapp W, Dörken B, Gilks WR, et al, eds. *Leukocyte Typing IV: White Cell Differentiation Antigens*. New York, NY: Oxford University Press; 1989:342-343.
27. Gallagher PF, Fazekas de St. Groth B, Miller JFAP. CD4 and CD8 molecules can physically associate with the same T-cell receptor. *Proc Natl Acad Sci, USA*. 1989;86:10044-10048.
28. Anderson P, Blue M-L, Morimoto C, Schlossman SF. Cross-linking of T3 (CD3) with T4 (CD4) enhances the proliferation of resting T lymphocytes. *J Immunol*. 1987;139:678-682.
29. Eichmann K, Jönsson J-I, Falk I, Emmrich F. Effective activation of resting mouse T lymphocytes by cross-linking submitogenic concentrations of the T cell antigen receptor with either Lyt-2 or L3T4. *Eur J Immunol*. 1987;17:643-650.
30. van Dongen JJM, Krissansen GW, Wolvers-Tettero ILM, et al. Cytoplasmic expression of the CD3 antigen as a diagnostic marker for immature T-cell malignancies. *Blood*. 1988;71:603-612.
31. Brenner MB, McClean J, Dialynas DP, et al. Identification of a putative second T cell receptor. *Nature*. 1986;322:145-149.
32. Clevers H, Alarcón B, Wileman T, Terhorst C. The T cell receptor/CD3 complex: a dynamic protein ensemble. *Annu Rev Immunol*. 1988;6:629-662.
33. Centers for Disease Control. Update: universal precautions for prevention of transmission of human immunodeficiency virus, hepatitis B virus, and other bloodborne pathogens in health-care settings. *MMWR*. 1988;37:377-388.
34. *Protection of Laboratory Workers from Infectious Disease Transmitted by Blood, Body Fluids, and Tissue: Tentative Guideline*. Villanova, PA: National Committee for Clinical Laboratory Standards; 1991. NCCLS document M29-T2.
35. Nicholson JK, Browning SW, Orloff SL, McDougal JS. Inactivation of HIV-infected H9 cells in whole blood preparations by lysing/fixing reagents used in flow cytometry. *J Immunol Methods*. 1993;160:215-218.
36. *Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture: Approved Guideline*. Villanova, PA: National Committee for Clinical Laboratory Standards; 1991. NCCLS document H3-A3.
37. *Clinical Applications of Flow Cytometry: Quality Assurance and Immunophenotyping of Peripheral Blood Lymphocytes; Tentative Guideline*. Villanova, PA: National Committee for Clinical Laboratory Standards; 1992. NCCLS document H42-T.
38. Giorgi JV. Lymphocyte subset measurements: significance in clinical medicine. In: Rose NR, Friedman H, Fahey JL, eds. *Manual of Clinical Laboratory Immunology*. 3rd ed. Washington, DC: American Society for Microbiology; 1986:236-246.
39. Jackson AL, Warner NL. Preparation, staining, and analysis by flow cytometry of peripheral blood leukocytes. In: Rose NR, Friedman H, Fahey JL, eds. *Manual of Clinical Laboratory Immunology*. 3rd ed. Washington, DC: American Society for Microbiology; 1986:226-235.

40. Prince HE, Hirji K, Waldbeser LS, Plaeger-Marshall S, Kleinman S, Lanier LL. Influence of racial background on the distribution of T-cell subsets and Leu 11-positive lymphocytes in healthy blood donors. *Diagn Immunol*. 1985;3(1):33-37.
41. Angadi CV. Lack of Leu-3a epitope on T-helper (CD4) lymphocytes. *J Clin Lab Anal*. 1990;4:193-195.
42. *How to Define and Determine Reference Intervals in the Clinical Laboratory: Approved Guideline*. Villanova, PA: National Committee for Clinical Laboratory Standards; 1992. NCCLS document C28-A.
43. Lanier LL, Le AM, Phillips JH, Warner NL, Babcock GF. Subpopulations of human natural killer cells defined by expression of the Leu-7 (HNK-1) and Leu-11 (NK-15) antigens. *J Immunol*. 1983;131:1789-1796.

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BD, логотип BD и другие товарные знаки являются собственностью компании Becton, Dickinson and Company. © 2010 BD